

How to use these Therapy Programs:

The **LazrPulsr 4X Mk III Home Therapy Support Laser** has multiple therapy programs built into the laser. Each program is designed to aid in the recovery of a specific health condition. Each therapy program consists of one or more preprogrammed frequencies that turn each laser beam on and off a specific number of times per second. You will not be able to see most of these pulses as most of them are much faster than the eyes can see.

Each different program is designed to target different tissues or stimulate specific healing actions while assisting in management of pain. A complete therapy protocol for any given condition consists of one to several different but related programs. These other recommended programs are listed to the side of each program.

Home Therapy Cards:

There are no Therapy Cards for the Home Therapy System. The homeopathic remedies for these programs are already imprinted into the laser itself.

Selecting a Therapy Program:

Each procedure is numbered at the top left corner of the image. Enter the number into the laser as follows.

1. Press the "Select" key and the program # is highlighted.
2. Enter the desired program number with the key pad.
3. Press the "Select" key again and the program is ready to operate.

Operating the program:

The Laser unit has been programmed to run an entire Therapy Program with one button push. Each program is broken down into Steps. Each Step is made up of one frequency program that will operate for a predetermined length of time then pause before beginning the next Step.

You may start by pressing the "Laser On/Off" key.

Pressing the "Laser On/Off" key will cause the laser to run the entire Therapy Program and pause for several seconds between each Step that requires a change in the laser's location on the body and then automatically continues the next Step.

You may start by pressing the "Time" key: (Highly Recommended)

Pressing the "Time" key, highly recommended, will cause the laser to run the entire Therapy Program but pause between each Step that requires the operator to change the laser's position on the body. The following Step number will be displayed while paused so the operator can locate the next area to be treated from this instruction guide. To continue with the next step press the "Laser Pulse" key. The operator can also pause the program at any time for up to 5 minutes by pressing the "Laser Pulse" key and continue by pressing the "Laser Pulse" key again.

You may start in the middle of a program:

If you desire to run one specific step of a program or start in the middle of a program that you began and did not have time to complete, follow these directions. Start the program with the "Time" key. Pause the program with the "Laser Pulse" key to pause the program. Press the Up or Down arrow keys to move to the desired step. Press the "Laser Pulse" to continue the program from that step.

You may run any program all in one location:

In some cases it may be much easier for the operator to just place the unit in one location to perform the entire procedure. If this is the case you will have to perform the procedure at least two times to accomplish the same therapeutic effect. The body location for this is indicated at the bottom left corner of each card.